

Dual-Source Filtration for Purity & Safety: A Research Framework for Medical-Grade Cosmeceutical Serums

Issue date: 2025-August-30

SM DERMA R&D Department

ARTICLE INFO	ABSTRACT
<p><i>Keywords:</i></p> <p>Exosome</p> <p>Stem Cell</p> <p>Medical-Grade</p> <p>Cosmeceutical</p>	<p>Extracellular vesicles (EVs) from human stem-cell cultures and plant sources show promise for dermal regeneration and antioxidant defense, but translation requires rigorous purification and safety controls.</p> <p>Objective: To outline and justify a dual-source workflow that isolates nanosized vesicles (50–150 nm) from human mesenchymal stem-cell conditioned media (hMSC-CM) and botanical extracts using ultrafiltration (UF) followed by size-exclusion chromatography (SEC), combined under sterile, endotoxin-free (human EVs) and pesticide-free (plant EVs) specifications. Methods/Significance: The platform aligns with MISEV2023 guidance on EV production, separation, characterization, and reporting, and leverages recent advances demonstrating SEC's advantages for EV purity and reproducibility. Outcomes (anticipated): Hybrid human–plant vesicle serums for post-procedure recovery and photoaging support, backed by preclinical potency assays and dermatologist-relevant endpoints.</p>

1. Introduction

EVs are nanoscale, lipid-bilayer vesicles that carry microRNAs, proteins, and lipids with roles in tissue repair and immune modulation. hMSC-EVs accelerate wound closure, enhance collagen remodeling, and reduce inflammation in skin models and early clinical translation, positioning them as cell-free candidates for regenerative dermatology. Plant-derived exosome-like

nanoparticles (PELNs) provide biocompatible carriers rich in antioxidant/anti-inflammatory phytochemicals and can aid skin delivery. Standardized isolation and safety testing are essential to avoid protein/lipoprotein co-isolation, residual endotoxin, or agricultural contaminants. [BioMed Central +3 PubMed Central +3 PubMed +3](#)

Hypothesis. A dual-source, orthogonal purification (UF→SEC) can yield high-purity hMSC-EVs and PELNs that are analytically defined, reproducible, and safe for medical-grade cosmeceutical applications. This design follows MISEV2023 recommendations emphasizing transparent reporting of production, separation, and functional assays. [ISeV Journals +1](#)

2. Rationale and Prior Evidence

2.1 Human hMSC-EVs for skin regeneration

Systematic and narrative reviews show hMSC-EVs promote angiogenesis, fibroblast proliferation, and remodeling with favorable safety compared to cell therapy; recent meta-analyses and umbrella reviews extend efficacy signals across wound repair contexts. [PubMed Central +2 Wiley Online Library +2](#)

2.2 Plant-derived vesicles (PELNs) as protective nanocarriers

Recent updates describe PELNs from citrus, ginger, grape, and green tea as stable, non-toxic vesicles that carry intrinsic antioxidants, enhance barrier function, and may facilitate cutaneous delivery of co-actives. [MDPI +1](#)

2.3 Why UF→SEC?

SEC increasingly replaces precipitation/ultracentrifugation for applications demanding higher EV purity and functional reproducibility, offering superior removal of protein contaminants and lipoproteins across biofluids. [ScienceDirect +2 Frontiers +2](#)

3. Methods (Proposed Workflow)

3.1 Source materials

Human: Serum-free hMSC-CM (mycoplasma-negative cell banks, early passages; chemically defined media).

Plant: Food-grade or pharmacopeial botanicals (e.g., *Camellia sinensis*, Citrus spp., *Zingiber officinale*), washed and cold-pressed/juiced under HACCP.

3.2 Purification sequence

1. Pre-clear: 0.45 μm → 0.22 μm filtration of CM/plant juice to remove cells/debris.
2. Tangential-flow UF: 100–300 kDa membranes to concentrate vesicles and remove <10 nm solutes.
3. SEC: Calibrated porous media to collect 50–150 nm EV fractions, minimizing co-isolated proteins/lipoproteins.
4. Sterile formulation: Isotonic buffer (pH 6.0–7.0), oxygen-barrier packaging; 2–8 °C storage (stability per ICH). Method reporting will follow MISEV2023 checklists (inputs, yields, fraction IDs). [ISeV Journals +1](#)

3.3 Characterization & release

- Identity/size: NTA/TRPS (mode 50–150 nm), cryo-TEM; human EV markers (CD9/CD63/CD81) and negative controls; plant vesicle markers per PELN literature.
- Purity: Particle:protein ratio; ApoB/lipoprotein depletion indices; RNA content.
- Potency (in vitro): Fibroblast scratch closure, keratinocyte migration, VEGF/IL-10 modulation, UV-ROS suppression in 3D skin models.
- **Safety:**
 - ✓ Endotoxin: LAL (chromogenic) with low endotoxin recovery (LER) controls and TLR4-blocking in bioassays, acknowledging masking risks documented for LPS in formulated matrices. [ScienceDirect +2 PubMed +2](#)
 - ✓ Pesticides (plant): Targeted LC–MS/MS panel vs. MRLs; require non-detect (nd) for release.

✓Bioburden/sterility:Per cosmetic-device context; mycoplasma-free certification for hMSC production.

✓Irritation/sensitization:HR IPT and in vitro cytokine panels.

4. Safety-by-Design and Quality System

- Endotoxin-free (human EVs):Controlled upstream culture, endotoxin-screened raw materials, validated cleaning; spike-and-recovery/LER investigations to prevent false negatives; JoVE protocols provide practical controls for endotoxin avoidance and evaluation. [PubMed +1](#)

- Pesticide-free (plant EVs):Certified supply chains; pre-harvest/COA verification; SEC helps reduce low-MW contaminants; LC–MS confirmation on final bulk. [MDPI](#)

- Documentation:Batch records include NTA counts, mode size, purity indices, endotoxin (EU/mL), pesticide panel, sterility/bioburden, and potency readouts—MISEV-aligned data package for each lot. [PubMed Central](#)

5. Experimental Plan

5.1 Preclinical

✱In vitro:

- Human dermal fibroblasts: scratch closure (%/24–48 h), COL1A1/COL3A1 qPCR, pro-collagen ELISA.

- Keratinocytes: migration, IL-1 α /IL-8 after UV-A; antioxidant endpoints in ROS probes.

- 3D skin equivalents: TEER recovery, histology (H&E, Masson), OCT microvasculature.

✱In vivo (non-clinical):

- Murine tape-stripping & UV-photodamage: TEWL, erythema, histologic collagen ratios.

- Safety: local tolerance; HRIPT (human) for finalized serum.

Justification:Cutaneous benefit signals for hMSC-EVs and plant vesicles are well documented across wound and photodamage models, enabling power calculations for effect sizes (e.g., wound closure and dermal density). [PubMed Central](#)

5.2 Early clinical (cosmetic endpoints)

Randomized, split-face adjunct after fractional laser (n \approx 30): downtime (erythema/edema days), profilometry (Ra/Rz), ultrasound dermal density, blinded photo-scoring, and non-invasive cytokine tape-strips. This mirrors EV-based cutaneous studies and adheres to cosmetic-claim boundaries. [ScienceDirect](#)

6. Application Concept: Medical-Grade Cosmeceutical Serum

✱Target use:Post-procedure recovery and photoaging support.

✱Formulation:Aqueous serum (pH 6.0–7.0) with HA/panthenol; EV integrity preserved (no surfactants that disrupt vesicle membranes).

✱Dose guidance:Start at $1\text{--}5 \times 10^9$ particles/mL (per NTA) adjusted by in-vitro EC50; apply in clinic (3–4 sessions, 2–4 weeks apart) + home-use for 4–8 weeks; finalize after Phase 0/feasibility data. (SEC-based workflows report robust particle recovery with improved purity, supporting dose precision.) [Frontiers](#)

7. Regulatory & Ethics

✱Nomenclature & reporting:Use “extracellular vesicles (EVs)” per MISEV2023; avoid definitive “exosome” biogenesis claims unless demonstrated. [ISeV Journals](#)

✱Claims & path:Cosmetic (non-therapeutic) language for over-the-counter serums; consider device/drug routes if therapeutic claims are

pursued (requires RCTs and regional regulatory submissions).

※Human-cell sourcing:IRB/ethics approvals, donor consent, traceability, mycoplasma-free certification.

※Plant sourcing:HACCP/GACP principles; sustainability documentation.

8. Discussion

Innovation.The dual-source UF→SEC strategy provides orthogonal purification that increases vesicle purity and removes confounders (free proteins, lipoproteins), directly addressing reproducibility concerns raised in EV research and enabling credible cosmetic translation. [PubMed Central](#)

Mechanistic complementarity.hMSC-EVs provide pro-regenerative signals (angiogenesis, matrix remodeling), while PELNs contribute antioxidant/anti-inflammatory support and may enhance delivery of co-actives—yielding a broader, more resilient efficacy profile. [PubMed Central +1](#)

Risk controls.Endotoxin masking (LER) can produce false “clean” results; validated LAL with LER challenges and orthogonal bioassays (e.g., TLR4 blockade) are mandatory. Agricultural inputs demand pesticide screening with LC–MS; release requires non-detect. [ScienceDirect +1](#)

Limitations.EV heterogeneity persists despite SEC; standardized reporting and potency assays are needed. Plant vesicle composition varies by species and harvest; strict specification and COAs mitigate variability. Long-term stability and scale-up economics require further study. [ISeV Journals](#)

9. Conclusion

A dual-source filtrationplatform—UF concentration followed by SEC purification—can generate high-purity, safety-verifiedhuman

and plant EV preparations suited for medical-grade cosmeceutical serums. By adhering to MISEV2023standards, implementing endotoxin-freeand pesticide-freerelease criteria, and validating function with skin-relevant bioassays and early clinical endpoints, this approach offers a credible path from bench to clinic-adjacent skincare.

Key references

1. Mesenchymal stem cell-derived extracellular vesicles in skin wound healing and scar repair by Ding J.Y. et al. – Int J Mol Sci.2023. Open access review of MSC-EVs for skin regeneration. PubMed Central
2. Extracellular vesicles derived from mesenchymal stem cells: Mechanisms in tissue repair and anti-fibrosis by Kou M. et al. – Cell Death & Disease.2022. Addresses EVs’ roles in tissue repair including skin injury. Nature
3. Plant-derived exosome-like nanovesicles in chronic wound healing by Wu W. et al. – Int J Nanomedicine.2024. Focuses on plant-derived vesicles (PELNs) for wound healing applications. PubMed Central
4. Plant-derived exosome-like nanoparticles: Emerging nanosystems in regenerative medicine – Int J Mol Sci.(2024). Comprehensive review of PELNs for tissue engineering and regeneration.